

ゲルろ過用アガロースゲルのメーカー比較

Agarose %	GE Biosciences Limited™	Cross-linked	BioRad	Cross-linked	ABT	Cross-linked	Recommended OR*
2%	Sepharose™ 2B	No	BioGel A-50 M	No	2% PLAIN Agarose Beads	No	100 – 50,000 kDa
2%	Sepharose™ CL-2B	Yes	None	N/A	2% CROSSLINKED Agarose Beads	Yes	-
4%	Sepharose™ 4B	No	BioGel A-15M	No	4% PLAIN Agarose Beads	No	40-15,000 kDa
4%	Sepharose™ CL-4B	Yes	None	N/A	4% CROSSLINKED Agarose Beads	Yes	-
4%	Sepharose™ 4B	No	BioGel A-5M	No	4% PLAIN Agarose Beads	No	40-15,000 kDa
6%	Sepharose™ CL-6B	Yes	None	N/A	6% CROSSLINKED Agarose Beads	Yes	-
6%	None	N/A	BioGel A-1.5M	No	6% PLAIN Agarose Beads	No	10-1,500 kDa
8%	None	N/A	None	N/A	8% CROSSLINKED Agarose Beads	Yes	-
8%	None	N/A	BioGel A-0.5M	No	8% PLAIN Agarose Beads	No	10-500kDa
10%	None	N/A	None	N/A	10% CROSSLINKED Agarose Beads	Yes	-

* Recommended operation ranges

発売元
プロテノバ株式会社
〒769-2604
香川県東かがわ市西村1488-1
TEL: 0879-49-0702/FAX: 0879-49-0702
E-mail: info@protenova.com