

## ゲルろ過用アガロースゲルのメーカー比較

Agarose %	GE Biosciences Limited™	Cross-linked	BioRad	Cross-linked	ABT	Cross-linked	Recommended OR*
2%	Sepharose™ 2B	No	BioGel A-50 M	No	2% PLAIN Agarose Beads	No	100 - 50,000 kDa
2%	Sepharose™ CL-2B	Yes	None	N/A	2% CROSSLINKED Agarose Beads	Yes	-
4%	Sepharose™ 4B	No	BioGel A-15 M	No	4% PLAIN Agarose Beads	No	40 - 15,000 kDa
4%	Sepharose™ CL-4B	Yes	None	N/A	4% CROSSLINKED Agarose Beads	Yes	-
4%	Sepharose™ 4B	No	BioGel A-5 M	No	4% PLAIN Agarose Beads	No	40 - 15,000 kDa
6%	Sepharose™ CL-6B	Yes	None	N/A	6% CROSSLINKED Agarose Beads	Yes	-
6%	None	N/A	BioGel A-1.5 M	No	6% PLAIN Agarose Beads	No	10 - 1,500 kDa
8%	None	N/A	None	N/A	8% CROSSLINKED Agarose Beads	Yes	-
8%	None	N/A	BioGel A-0.5 M	No	8% PLAIN Agarose Beads	No	10 - 500 kDa
10%	None	N/A	None	N/A	10% CROSSLINKED Agarose Beads	Yes	-

\*Recommended operating ranges