



AFFINITY HIS-TAG PURIFICATION

Chelating Rapid Run™ Agarose Beads

PRODUCT	Metal Free Rapid Run™	Nickel Rapid Run™	Cobalt Rapid Run™
CAT.	6RR-QH-X	6NiRR-X	6CoRR-X
Bead geometry & size	Spherical, Standard: ~ 90 – 150 µm		
Exclusion limit	~4x10 ⁶		
Crosslinked	Highly Crosslinked		
Agarose%	6%		
Binding / Loading capacity (µmol Me ²⁺ / ml gel)	~ 20		
Ligand	Iminodiacetic Acid		
Antimicrobial Agent	20% ethanol		
Storage temperature	2 – 30°C		

X: Product quantity 25, 100 and 500 ml.
For laboratory use only. Not for use in diagnostic or therapeutic procedures.

ABT-Chelates-RR_TS-2022A