

AFFINITY COUPLING

Aminoethyl Agarose Resin

PRODUCT	Low-Density Aminoethyl 6BCL	Very Low-Density Aminoethyl 4BCL	High-Density Aminoethyl 6BCL	High-Density Aminoethyl 4BCL
Cat. No.	6BCL-AL0-X	4BCL-AVL4-X	6BCL-AM3-X	4BCL-AH1-X
Bead Geometry & Size	Spherical, Standard ~ 50 – 150 µm			
Crosslinked	Yes			
Matrix Active Groups	Amino Groups			
Agarose %	6%	4%	6%	4%
Activation Degree (µmol Diaminoethyl/MI Gel)	15 – 25	3 - 6	40 – 60	40 - 60
Antimicrobial Agent	20% Ethanol			
Storage Conditions	2 – 8 °C			

X: Product quantity. 25 ml or 100 ml.
For laboratory use only. Not for use in diagnostic or therapeutic procedures.

ABT-Aminoethyl-TS-2022A